



Low Linting



100% Cotton



X-Ray Detectable

ThreadX Packing Gauze is constructed from 100% cotton. Designed for open surgical procedures with an x-ray detectable thread heat bonded to the gauze.

ABSORBENT

Absorbs fluid and protects organ and tissue injury and helps control bleeding

X-RAY DETECTABLE

X-ray detectable thread for patient safety

LOW LINTING

Aids in infection control and minimises the risk of contamination



Suitable for:

Open surgical procedures

Deep wounds

Undermining and tunnels

Reorder	Ply	Size	Thread Colour	Sterility	Pcs/Unit	Units/Ctn
XRP001	4ply	1.25cm x 1m	Blue	Sterile	100	1
XRP006	4ply	2.5cm x 1m	Blue	Sterile	100	1
XRP007	4ply	2.5cm x 2m	Blue	Sterile	100	1
XRP011	4ply	5cm x 1m	Blue	Sterile	100	1
XRP012	4ply	5cm x 2m	Blue	Sterile	100	1
XRP016	4ply	10cm x 1m	Blue	Sterile	100	1
XRP017	4ply	10cm x 2m	Blue	Sterile	100	1
XRV003	8ply	10cm x 3m	Black	Sterile	100	1