

Leukomed®

LEUKOMED®
T PLUS

RELIABLE PROTECTION
AGAINST BACTERIA
AND VIRUSES

THERAPIES. HAND IN HAND.

BSN medical

LEUKOMED® T PLUS TRANSPARENT FILM DRESSING WITH ABSORBENT PAD

This water, bacteria and virus proof film dressing combines the ultra-thin and flexibility of the **Leukomed®** films with a low adherent wound contact layer. The low profile film and the strong adhesion allows the dressing to stay in place for longer enabling less dressing changes and decreased risk of wound contamination ensuring an optimal environment for wound healing.

EASY TO APPLY WITH UNIQUE RED TAB TECHNOLOGY

Enables easy and precise application for the right dressing placement

WATER, BACTERIA AND VIRUS PROOF

Protects reliably against water, bacteria and viruses to reduce infection risk. The waterproof barrier allows the patient to continue on with activities of daily living without the need for dressing changes, saving on time and additional product costs



COMFORTABLE AND CONFORMABLE TRANSPARENT FILM

Low profile film enables high wearing comfort for the patient

ABSORBENT AND BREATHABLE

Absorbent wound pad enables atraumatic removal and reduced dressing changes

STAYS IN PLACE

Rounded edges reduces the occurrence of roll back and snagging, allowing for a proper and secure hold

ORDERING INFORMATION

Ref-No.	Size	Pad Size	Items per Unit
Leukomed® T Plus			
72382-00	5 x 7.2 cm	2.8 x 4 cm	50
72382-01	8 x 10 cm	4 x 6 cm	50
72382-02	8 x 15 cm	4 x 11 cm	50
72382-03	10 x 25 cm	5 x 20 cm	50
72382-05	10 x 35 cm	5 x 30 cm	50

Indications

Leukomed® T plus – protection plus absorption

- For surgical incisions
- To cover and protect superficial cuts, lacerations and minor burns

BSN medical Australia & New Zealand

L6, 211 Wellington Road,
Mulgrave, Victoria, Australia 3170

www.bsnmedical.com.au
www.bsnmedical.co.nz

BSNM0050

04/2016

Customer Service Australia
T 1300 276 633 F 1300 998 830

Customer Service New Zealand
T 0508 276 111 F 0508 998 830

BSN medical